
Elettrodi Laringei Adesivi - Istruzioni per l'uso.....	1
Laryngeal Adhesive Electrodes - Instruction for use.....	9
Adhesive Laryngalelektrode - Gebrauchsanweisung.....	17
Électrodes Laryngées adhesives - Instructions.....	25
Electrodos Laríngeos Adhesivos - Instrucciones de uso.....	33
Eletrodos de laringe - Instrução de uso.....	41
喉返神经电极 - 使用说明书.....	49

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Laryngeal Adhesive Electrodes

Instruction for use

Laryngeal Electrodes for intraoperative monitoring.

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Manufacturer: Spes Medica S.p.A., Via Buccari 21, 16153 Genoa, Italy

All rights are reserved for this user manual, in particular the right to reproduction and distribution as well as translation. No part of this user manual may be reproduced in any form (e.g. photocopy, microfilm or other procedure) or processed, copied or distributed using electronic system without the prior written consent of Spes Medica.

CONTENTS

1. Intended Use
2. Handling
 - 2.1 Laryngeal Electrode Packaging
 - 2.2 Instructions
3. Waste disposal
4. Warnings
5. Symbols
6. Technical specifications
 - Connectors and cables
7. Manufacturer

1. Intended use

Spes Medica Laryngeal Adhesive Electrodes are electrodes to be used for the intraoperative monitoring of the larynx.

This product is a class II passive medical device to be used for intraoperative monitoring during surgical operation such as: Cervical Anterior Fusion, Carotid Endarterectomy, Thyroidectomy, Cranial Base procedures, Craniotomy surgery.

Note: The Laryngeal Electrode is a medical device to be used by specialized personnel.

2. Handling

The reliability of the muscular signal recorded through the Laryngeal Electrode is closely related both to the right electrode positioning on the endotracheal tube and to the accuracy of the endotracheal tube position in the patient's larynx.

2.1 Laryngeal electrode Packaging

Each Adhesive Laryngeal Electrode is provided in a sterile pouch.

Each box includes 10 pouches of Laryngeal Electrodes.

If the packaging is damaged, dispose of the electrode without use it.

2.2 Instructions for Use

① Choose the appropriate electrode size: the correct size depends on the endotracheal tube size on which it'll be fixed. Below the table with sizes indication:




Electrode code	Tube internal diameter	Channels	Size
LSE3446MCP075	2.0 – 5.5 mm	1	Pediatric*
LSE3446MCA075	6.0 – 7.0 mm	1	Adult
LSE3463DCA200	6.0 – 7.5 mm	2	Adult
LSE4263DCA200	8.0 – 9.5 mm	2	Adult
LSE4246MCA075	7.5 – 10.0 mm	1	Adult




**If the pediatric laryngeal electrode is chosen, trim the electrode along the numbered line on the electrode (i.e. if a 3 mm endotracheal tube is used, trim the electrode along the line marked with number 3).*

Make sure that there are no sharp angles after trimming that could hurt the patient.

② Checked both the integrity of the pouch before opening and the integrity of the laryngeal electrode once opened the pouch.

- ③ Apply the Laryngeal Electrode following the instruction in the following tables:

Laryngeal Electrode – SINGLE CHANNEL		
STEP 1	STEP 2	STEP 3
		
<p>Remove the protective film of the adhesive part of the electrode and apply the electrode on the endotracheal tube as shown in the figure. Place the electrode about 1 cm away from the cuff overlapping the red line to the endotracheal tube guideline.</p>	<p>Firstly, wrap on the tube the widest flap of the electrode (part “1” in the figure). Make sure the electrode is well-adhered without any creases.</p>	<p>Then wrap on the tube the smallest flap of the electrode (part “2” in the figure) Make sure the electrode is well-adhered without any creases.</p>

Laryngeal Electrode – DUAL CHANNEL		
STEP 1	STEP 1	STEP 1
		
<p>Remove the protective film of the adhesive part of the electrode and apply the electrode on the endotracheal tube as shown in the figure. Place the electrode about 1 cm away from the cuff overlapping the red line to the endotracheal tube guideline.</p>	<p>Firstly, wrap on the tube the widest flap of the electrode (part “1” in the figure). Make sure the electrode is well-adhered without any creases.</p>	<p>Then wrap on the tube the smallest flap of the electrode (part “2” in the figure) Make sure the electrode is well-adhered without any creases.</p>

For *pediatric size* laryngeal electrodes: firstly, wrap the smallest (non-recording part) flap of the electrode, then wrap the widest. Make sure the electrode is well-adhered without any creases.

④ A small amount of lubricant can be applied on the electrode avoiding the conductive parts of the electrode.

- ⑤ Finally proceed with the standard intubation procedure, making sure that the recording areas of the electrode are in contact with the vocal cords (Figure 1).


Figura 1

- ⑥ Connect the red and blue Touch Proof connectors of the laryngeal electrode to the ports + and - of EMG recording equipment (Figure 2). Connect the ground electrode connector (not included) to the appropriate port of EMG recording equipment (Figure 2, green cable).


Figure 2

- ⑦ Follow the next steps as described in the user manual of EMG equipment used.

3. Waste disposal

After use, Spes Medica Laryngeal Electrodes must be disposed of in accordance with regulations in force for medical devices.

4. Warnings

The Adhesive Laryngeal Electrode is single use only not designed for the re-use. Attempts to clean and re-use this single use device expose patients and operators to a risk of cross contamination and malfunction

Do not resterilize.

Avoid long-term using of paralyzing agents that could reduce or even delate the EMG response to the stimulation.

Avoid using of lubricant directly on the conductive part of the electrode.

Make sure the electrode is well-attached to the endotracheal tube and there are no sharp edges which could harm the patient.

Avoid excessive tractions to the cables.

To avoid contamination of the device and therefore of the patient, the user must wear disposable gloves

Only to be used by the clinical specialist.

It is necessary for the user and/or patient to report any serious incident occurring, related to the device, to the manufacturer and the competent authority of the member state where the user and/or patient is established.

The environment of use is a hospital, specifically an operating room.


The device is intended for both adult and pediatric patient populations. The device does not have limitations of use related to age (from birth onward), gender, or weight. The clinical specialist will identify the correct size for the clinical need.

There are no contraindications or warnings for pregnant or lactating women.

5. Symbols and explanations

Symbols	Meaning
	Complies with the Regulation (EU) 2017/745
	Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices
	注册证号码: 国械注进 20182071758
	Swiss authorized representative
	Medical device
	Unique device identifier
	Catalogue number
	Lot number
	Consult instruction for use
	Warning
	Single Use Only
	Do not resterilize
	Expiration date

Ethylene Oxide Sterilization with single sterile barrier system
with protective packaging outside

Ethylene Oxide Sterilization with single sterile barrier system

Do not use if package is damaged

Manufacturer

All Spes Medica Laryngeal Electrodes have UDI (Unique Device Identification) barcode GS1-128.

6. Technical specifications

Connectors and cables

The Spes Medica Laryngeal Electrodes are equipped with PVC insulated connection cables and a 1.5 mm Touch Proof connector.

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Spes Medica S.p.A. Tel: +39 (0)10 390343