Disposable Insulated Needle Electrode

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Intended use

The Needle Electrodes are intended for use for the recording of biopotential signals with recording, monitoring and stimulation equipment for neurophysiology, examples include: Electromyography (EMG), Electroencephalography (EEG), intraoperative monitoring (IOM) and evoked potentials (EP)




Direction for use

- Do not use the selected needle if the pouch is damaged or if the product has expired.
- Open the pouch and remove needle from storage pouch;
- Connect the needle to the equipment in use.
- Remove the needle sheath;
- Introduce the needle into the tissue to be tested. Insert the needle up to hub;
- If the needle bends before, during or after insertion, do not straighten or re-insert.
- Seal and replace with a new one.
- After the use of the needle according to standard protocols or alternatively put the needle in approved biohazard sharps containers.

Warnings

- The resulted needle electrode is disposable.
- This product is not designed for re-use. Attempts to clean and re-use this single-use device expose patients and operators to a risk of cross contamination.
- Do not re-sterilize.

- For devices marked with  compatible with magnetic resonance imaging (MRI) is given as a range of applicability of 1.5 or 3.0 Tesla.
- To avoid contamination of the device and therefore of the patient, the user must wear disposable gloves
- Only to be used by the clinical specialist.
- It is necessary for the user and/or patient to report any serious incident occurring, related to the device, to the manufacturer and the competent authority of the member state where the user and/or patient is established.
- The environment of use is a hospital, specifically an operating room.
- The device is intended for both adult and pediatric patient populations. The device does not have limitations of use related to age (from both onward), gender, or weight. The clinical specialist will identify the correct size for the clinical need.
- There are no contraindications or warnings for pregnant or lactating women.
- Do not bend the needle and do not use instruments to modify its shape.
- Do not pull the electrode by the cable under any circumstances. Always grasp the handle for insertion and removal.




	Complies with the Regulation (EU) 2017/745
Rx Only	Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices
	Swiss authorized representative
MD	Medical device
UDI	Unique device identifier
REF	Catalogue number
LOT	Lot number
	Single Use Only
	Do not reuse if package is damaged
	Consult instruction for use
	Do not use if package is damaged
	Warnings and Cautions associated with this device.
	Expire by
	Ethylene Oxide Sterilization with single sterile barrier system
	Ethylene Oxide Sterilization with single sterile barrier system with protective packaging outside
	Manufacturer
	MRI Conditional

Disposable Insulated Needle Electrode Instructions for use

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