
Elettrodi Laringei con tubo endotracheale integrato	1
Istruzioni per l'uso	
Laryngeal Electrodes with integrated E.T. tube	8
Instruction for use	
Laryngal Elektrode mit integriertem E.T. Tubus	15
Gebrauchsanweisung	
Électrodes Laryngées avec sonde d'intubation intégrée	22
Instructions	
Electrodos laríngeos con tubo endotraqueal integrado	29
Instrucción de utilizo	
Eletrodo De Laringe Com Tubo Enotraquel Integrado	36
Instruções de uso	

ITALIANO

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Laryngeal Electrodes with integrated E.T. tube

Instruction for use

Laryngeal Electrodes for intraoperative monitoring.

Manufacturer: Spes Medica S.p.A., Via Buccari 21, 16153 Genoa, Italy

All rights are reserved for this user manual, in particular the right to reproduction and distribution as well as translation. No part of this user manual may be reproduced in any form (e.g. photocopy, microfilm or other procedure) or processed, copied or distributed using electronic system without the prior written consent of Spes Medica.

CONTENTS

1. Intended Use
2. Handling
 - 2.1 Laryngeal Electrode Packaging
 - 2.2 Instructions
3. Waste disposal
4. Warnings
5. Symbols
6. Technical specifications
 - Connectors and cables
7. Manufacturer

1. Intended use

Spes Medica Laryngeal Electrodes are surface electrodes to be used for the intraoperative monitoring of the larynx.

This product is a class II passive medical device to be used for intraoperative monitoring during surgical operation such as: Cervical Anterior Fusion, Carotid Endarterectomy, Thyroidectomy, Cranial Base procedures, Craniotomy surgery.

The Laryngeal Electrode is already placed on the armed or not armed endotracheal tube.

Note: The Laryngeal Electrode is a medical device to be used by specialized personnel.

2. Handling

The reliability of the muscular signal recorded through the Laryngeal Electrode is closely related to the accuracy of the endotracheal tube positioning.

2.1 Laryngeal electrode Packaging

Each Laryngeal Electrodes is provided in a sterile pouch.

Each box includes 5 pouches of Laryngeal Electrodes.

If the packaging is damaged dispose of the Laryngeal Electrode without use it.

2.2 Instructions for Use

① Choose the appropriate electrode size: the Laryngeal Electrode with integrated E.T. tube has the correct size matches with the size of endotracheal tube to be used. The dimension of the tube, internal diameter, is shown on the product label (unit of measure mm).

② Check both the integrity of the pouch before opening, the integrity of the laryngeal electrode after the pouch opening and the shelf life.

③ Open the pouch from the Chevron side.

④ In case of need for lubrication, apply the lubricant on the endotracheal tube avoiding the recording part of the electrode.

⑤ Position correctly the patient: the head must be in hyperextension with the nape as close as possible to the cervical vertebrae, reaching the so-called “sniffing” position.

⑥ Insert the laryngoscope into the mouth.

⑦ Insert the endotracheal tube with integrated laryngeal electrode, if necessary with the stylet placed in the tube, pushing it through the vocal cords spread on the laryngoscope guide. Making sure that the recording areas of the electrode are in contact with the vocal cords (Fig. 1).



Fig. 1

- ⑧ If used, remove the stylet.
- ⑨ Remove the laryngoscope taking care to hold the position of the endotracheal tube with integrated laryngeal electrode.
- ⑩ Inflate the cuff with enough ml of air using a syringe without tip to allow the tube fixing to the trachea.
BE CAREFUL: excessive inflation of the cuff could lead to the ischemia of the tracheal mucosa reducing its bloody supply and damaging it.
- ⑪ Check the correct positioning of the electrode in trachea by manual ventilation.
- ⑫ Make sure no air comes out form nose or mouth.
BE CAREFUL: air flowing from nose or mouth indicates that the cuff is deflated or damage.
- ⑬ Fasten the endotracheal tube with two pieces of plaster wrapping them around the tube and then fixing to the upper part of the lip and to the cheeks.

(Fig. 2)

- ⑭ Finally start the mechanical ventilation.
- ⑮ Connect the red and blue Touch Proof connectors of the laryngeal electrode to the ports + and - of EMG recording equipment (Fig. 3).
Connect the ground electrode connector (not included) to the appropriate port of EMG recording equipment (Fig. 3, green cable).


Fig. 3

- ⑯ Follow the next steps as described in the user manual of EMG equipment used.

3. Waste disposal

Spes Medica Laryngeal Electrodes must be disposed of in accordance with regulations in force for medical devices.

4. Warnings

The Laryngeal Electrode with E.T. integrated tube is single use only not designed for the re-use. Attempts to clean and re-use this single use

device expose patients and operators to a risk of cross contamination and malfunction.

Do not resterilize.

Avoid using of paralyzing agents that could reduce or even delate the EMG response to the stimulation.

Avoid using of lubricant directly on the conductive part of the electrode.

Do not over inflate the balloon. Over inflation may lead to the balloon breaking causing damages to the patient's tracheal and/or the tube blockage.

Make sure the balloon is totally flat before removing it from the tracheal.

Do not perform MRI exams on patient intubated with laryngeal electrodes. The effect of MRI on the reinforcing coiled wire of the tube is unknown.



During the monitoring do not rotate the patient's neck and head, this could lead to the wrong electrode placement and to an altered response.

To avoid contamination of the device and therefore of the patient, the user must wear disposable gloves

Avoid the Laryngeal Electrodes with an integrated endotracheal tube in the direction of the nose, going to make an "S" bend, as shown in the figure. This can cause malfunction and reduce the performance of the electrode

Only to be used by the clinical specialist.

	It is necessary for the user and/or patient to report any serious incident occurring, related to the device, to the manufacturer and the competent authority of the member state where the user and/or patient is established.
	The environment of use is a hospital, specifically an operating room.
	The device is intended for both adult and pediatric patient populations. The device does not have limitations of use related to age (from birth onward), gender, or weight. The clinical specialist will identify the correct size for the clinical need.
	There are no contraindications or warnings for pregnant or lactating women.

5. Symbols and explanations

Symbols	Meaning
	Complies with the Regulation (EU) 2017/745
	Swiss authorized representative
	Medical device
	Unique device identifier
	Catalogue number
	Lot number
	Consult instruction for use
	Warning

	Single Use Only
	Do not resterilize
	Expire date
	Ethylene Oxide Sterilization with single sterile barrier system with protective packaging outside
	Ethylene Oxide Sterilization with single sterile barrier system
	Do not use if package is damaged
	Manufacturer

All Spes Medica Laryngeal Electrodes have UDI (Unique Device Identification) barcode GS1-128.

6. Technical specifications

Connectors and cables

The Spes Medica Laryngeal Electrodes are equipped with PVC insulated connection cables and a 1.5 mm Touch Proof connector.

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