




Bridge Electrodes / EEG

Bridge Electrodes

- Gold screw guarantees high signal quality.
- Screw system suitable for correct positioning of the electrode on the skin.
- Universal connector suitable for 2 mm spring pin and crocodile clip.
- Bridge made of ABS plastic material.



Pack of 10 connection cables are available in the following colors:



| Code | Electrode Material | Size | Bridge Color | Pack |
|----------|--------------------|-----------|--------------|--------|
| BEAG00A2 | AgAgCl | Adult | Red | 10 pcs |
| BEAG00A7 | AgAgCl | Adult | Blue | 10 pcs |
| BEAG00P4 | AgAgCl | Pediatric | Yellow | 10 pcs |
| BEAG00N1 | AgAgCl | Neonatal | No Bridge | 10 pcs |
| BEG000A2 | AgAgCl Sinterised | Adult | Red | 10 pcs |
| BEG000A7 | AgAgCl Sinterised | Adult | Blue | 10 pcs |
| BEG000P4 | AgAgCl Sinterised | Pediatric | Yellow | 10 pcs |

Connection Cables

- Cable is Tinsel-coated guaranteeing high flexibility.
- Silver conductors provide maximum mechanical resistance and excellent electrical conductivity.

| Code | Electrode Side | Headbox Side | Pack |
|-------------|---------------------|-----------------|--------|
| BEC1226C1GS | Insulated Crocodile | Ø 1,5 mm | 10 pcs |
| BEC1223C1GS | Insulated Crocodile | Spring 2 mm pin | 10 pcs |
| BEC1221C1GS | Insulated Crocodile | 2 mm Pin | 10 pcs |
| BEC1226C1PS | Mini Croco-Clip | Ø 1,5 mm | 10 pcs |
| BEC122326S | Spring 2 mm pin | Ø 1,5 mm | 10 pcs |
| BEC122323S | Spring 2 mm pin | Spring 2 mm pin | 10 pcs |
| BEC122631S | Spring 3 mm pin | Ø 1,5 mm | 10 pcs |
| BEC122641S | Spring 4 mm pin | Ø 1,5 mm | 10 pcs |

Cable Length 120 cm

Assembling Kit for Bridge Electrodes

| Code | Description | Pack |
|----------|---|---------|
| BCSET000 | Cotton pads, rubber rings and covering tool | 100 pcs |

Bridge Electrode Bowl

| Code | Description | Pack |
|-----------|-----------------------|------|
| EHOS00000 | Bridge Electrode Bowl | 1 pc |