




PSG / Sensors


Reusable Airflow Thermal Sensor

- Overmoulded in soft plastic probes for oronasal breathing recording.
- Made of flexible materials developed to ensure optimized patient comfort.

Code	Model	Size	Pack
AIRFA00300	Airflow 3 ways	Adult	1 pc
AIRFP00300	Airflow 3 ways	Pediatric	1 pc
AIRFN00300	Airflow 3 ways	Neonatal	1 pc
 Touch Proof Ø 1,5 mm DIN 42802 Cable Length 200 cm Available in 1 way, 2 ways or divided ways			



Disposable Airflow Thermal Sensor


Code	Model	Size	Pack
AIRFA0030D	Airflow 3 ways	Adult	3 pcs
 Touch Proof Ø 1,5 mm DIN 42802 Cable Length 200 cm			



Upgrade +

Reusable Respiratory Effort Sensor


- Band is made with elastic and comfortable materials for increased patient comfort.
- Piezoelectric sensor.

Code	Model	Size	Pack
RESPA00000	Abdomen/Chest	Adult	1 pc
RESPP00000	Abdomen/Chest	Pediatric	1 pc
 Touch Proof Ø 1,5 mm DIN 42802 Cable Length 200 cm			



Disposable Respiratory Effort Sensor

- No need chest belt.
- Comfortable for children.

Code	Model	Size	Pack
RESPA000200	Abdomen/Chest	Adult	5 pcs
RESPP000200	Abdomen/Chest	Pediatric	5 pcs
 Touch Proof Ø 1,5 mm DIN 42802 Cable Length 200 cm			


Code	Description	Cable Length	Pack
BMDS00260R	Connection Cable for Disposable Respiratory Effort Sensor	60 cm	1 pc






Tremor Sensor


- Suitable for tremor analysis.
- The device is reusable.


Code	Recorded Signal	Cable Length	Pack
TREM000000	Tremor	200 cm	1 pc
 Touch Proof Ø 1,5 mm DIN 42802			



Limb Movement Sensor


- Excellent patient comfort and long-lasting product over time.
- The device is reusable.

Code	Recorded Signal	Cable Length	Pack
LIMB000000	Limb movement	200 cm	1 pc
 Touch Proof Ø 1,5 mm DIN 42802			




Reusable Snoring Sensor

- Converts snoring noises into an analogue signal.
- Ensures high quality recording in noisy environments.

Code	Recorded Signal	Cable Length	Pack
SNOR000000	Snoring	200 cm	1 pc
 Touch Proof Ø 1,5 mm DIN 42802			



Disposable Snoring Sensor

Code	Recorded Signal	Cable Length	Pack
SNOR000200	Snoring	200 cm	5 pcs
 Touch Proof Ø 1,5 mm DIN 42802			

Code	Description	Cable Length	Pack
BMDS00260S	Connection Cable for Snoring Sensor	60 cm	1 pc