




PSG / Bands, Sensor and Airflow Cannulas

Inductive Band


- Made with comfortable materials and equipped with an adjustable buckle for increased patient comfort.
- Easily washable.
- The device is reusable.

Code	Model	Connectors	Pack
9010	Inductive Band	Snap	1 pc
9101	Connection Cable	 Ø 1,5 mm	1 pc



Position Sensor

- The "Position Sensor" kit is suitable for PSG recordings in the following positions: prone, supine, left side, right side, standing and sitting.
- Includes an adjustable elastic band.

Code	Sensor	Pack
SSP00026AC	Position Sensor AC	1 pc
SSP00026AC-kit	Position Sensor AC	1 pc
SSP00026DC	Position Sensor DC	1 pc
SSP00026DC-kit	Position Sensor DC	1 pc
 Touch Proof Ø 1,5 mm DIN 42802		
Cable Length 200 cm		



Thermocan

- Thermocan is designed with integrated thermistor that allows a high quality recording of thermal flow and pressure monitoring.
- The device is disposable and complies with current AASM standards for nasal and oro/nasal thermal flow monitoring.

Code	Size	Cable Length	Pack
15807	Adult	210 cm	5 pcs
15807-2	Adult	60 cm	5 pcs
Code	Description	Cable Length	Pack
1475	Connection Cable	200 cm	1 pc



Nasal Cannula

Code	Size	N° Way	Cable Length	Filter	Pack
AIRFA2120D	Adult	2	210 cm	No	1 pc
AIRFA0620D	Adult	2	60 cm	No	1 pc
AIRFA212FD	Adult	2	210 cm	Yes	1 pc
AIRFA062FD	Adult	2	60 cm	Yes	1 pc

