



ECG

CLAMP ELECTRODE



UMDNS

11439

Description

Spes Medica has many years of experience in the production of reusable ECG electrodes.

In the various stages of production are employed high quality materials to ensure the best recording performance and durability of the product.

The electrodes are characterized by a universal connector and are sold in sets of 4 colors.

- A sensor 20% larger increases the recording capability and gives to the electrode the ability to record even in a bad placement;
- No plastic nut, only strong metal screwed nut to avoid any breaking point.

Materials

Sensor	AgAgCl, Ni
Clamp	Pa, Pom - C
Packaging	LDPE

Safety

Biocompatibility <i>(UNI EN ISO 10993)</i>	Yes
Latex	No
Phthalates	No

Specifications

Generic code	CLXXXXXX
Electrode dimensions	22x39 mm
Size	Adult, Pediatric
Electrode colour	● ● ● ●
Connector	Universal connection

Environment

PVC-Free	✓
RoHS II	✓
REACH	✓

Different configurations available on request