



ECG

GOLD ELECTRODES



UMDNS

11439

Description

They are made of high quality materials.

The recording surface is golded in order to have the best signal available.

The connector is overmoulded with the lead wire and it guarantees the maximum product duration.

Materials

Sensor	AgAgCl
Lead wire	Cu/Sn with PVC jacket
Connector	Cu/Au with PVC jacket
Packaging	PE

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	218XXXXXXX
Electrode diameter	12 - 24 mm
Connector	4 mm socket, Touch proof DIN 42802

Environment

RoHS II	✓
REACH	✓

Different configurations available on request