



EEG

MEDCAP Bonnets à Électrodes Précablées



UMDNS

17554

Description

Les principales caractéristiques de ces bonnets sont la fiabilité éprouvée depuis de nombreuses années avec tous les types d'équipements EEG, la possibilité de réparer rapidement et facilement ainsi que la grande flexibilité de ce type de produit pour le laboratoire EEG.

- Bonnet élastique;
- Code couleur pour identifier les tailles;
- Nettoyage et remplacement facile des électrodes.

Matériaux

Électrode	Étain (sur demande AgAgCl, AgAgCl sintérisé ou Doré)
Support	LDPE
Cordon	PVC
Bonnet	PU
Packaging	PP

Sécurité

Biocompatibilité (UNI EN ISO 10993)	Oui
Latex	Non
Nickel	Non
Phthalates	Non

Spécifications

Code Générique	CAMXXXXXXX
Nombre d'Électrodes	De 9 à 256
Connecteur	Disponible pour tous types de connexion

Environnement

RoHS II	✓
REACH	✓

Différentes configurations disponibles sur demande

Adult

head circumference about cm						
62 58	60 56	58 54	56 52	54 50	52 48	50 46
LARGE Cod. LA	LARGE/ MEDIUM Cod. LMA	MEDIUM Cod. MA	SMALL/ MEDIUM Cod. MSA	SMALL Cod. SA	EXTRA SMALL/ SMALL Cod. SEA	EXTRA SMALL Cod. SA
blue	blue	red	red	yellow	yellow green	green

Infant

head circumference about cm				
46 42	44 40	42 38	40 36	38 34
LARGE Cod. LI	LARGE/ MEDIUM Cod. LMI	MEDIUM Cod. MI	MEDIUM /SMALL Cod. MSI	SMALL Cod. SI
light blue	light blue	pink	pink	brown

Neonatal

head circumference about cm		
34 30	30 26	26 22
LARGE Cod. LN	MEDIUM Cod. MN	SMALL Cod. SN
yellow	blue	red