



SAC2



Description

SAC2 is an adhesive cream to be used with EEG Surface Electrodes. Provides resistance, cohesion and adhesiveness during use even if long term exams are made. SAC2 can be easily and quickly removed with water.

- Quick-drying cream
- No electrode damages
- Water soluble
- No residues
- No adverse smell
- External use only

UMDNS 11425

Composition

Chemical and physical properties

Specifications

Environment

Ingredients	Water, Talc, Celite, Glycerol, CarboxyMethylCellulose, Sodium chloride, Potassium chloride, Phenoxyethanol, Ethylhexylglycerin
Tube material	Aluminum
pH	8.2 ± 0.1
Generic code	SACX-XX
Packaging	10 tubes 100g each
RoHS II	✓
REACH	✓