




ADAPTER CABLE



Description

The adapter can stay connected to the headbox saving the headbox connector and the cap can be connected and disconnected using just the adapter.
Different configurations and connectors for the most of the EEG equipment's available on request.

UMDNS 17554

Materials

Lead wire	Cu/Ag with PVC jacket
Connector	CuZn/Au with PE jacket
Packaging	PE

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	CAPCABXXXXXX
Contacts number	From 10 to 36
Cable	Flat Twisted
Connector	Touch Proof 1,5 mm DIN 42802 D-SUB

Environment

RoHS II	✓
REACH	✓