




ADHESIVE ELECTRODES WITH REDUCED RECORDING AREA



Description

Spes Medica adhesive electrodes have a special shape and a specific size studied to allow a perfect registration even for the smaller muscles.

The electrode is made by biocompatible soft materials that give the best comfort to the patient.

- AgAgCl electrode for the best conduction;
- Easy snap connection;
- No gel on the electrode

UMDNS 11440

Materials

Safety

Specifications

Environment

Electrode	AgAgCl
Liner	PET
Connector	Stainless steel CuZn/Ni with PE jacket
Cable	Cu/Sn with PVC jacket
Packaging	AL/PE

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Generic code	DENIDXXXXX DAEDSXXXX
Recording area diameter	4 – 6 mm
Cable length	8 cm No cable
Connector	0,7 mm Touch Proof Pin Snap

RoHS II	✓
REACH	✓