




ADHESIVE ELECTRODES

DENIL



Description

Spes Medica adhesive electrodes have a special shape and a specific size studied to allow a perfect registration also on smaller muscles.

The innovative Hydrogel solid gel is developed especially for this application and gives to the electrode stable adhesion on the skin. It also can be repositioned several times.

- High quality recordings due to low impedance;
- Indicated for long-term testing;
- The electrode is made by biocompatible soft materials that give the best comfort to the patient;
- Available different shapes and sizes to allow a perfect registration also on smaller muscles.

UMDNS 11440

Materials

Safety

Specifications

Environment

Electrode	AgAgCl
Material in contact with the skin	Hydrogel
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Ni with PE jacket
Liner	PET
Packaging	AL/PE
Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No
Generic code	DENILXXXXX
Dimensions	20x25 mm
Cable length	From 8 to 250 cm
Cable colour	○ ● ● ●
Connector	Bipolar 0,7 mm Touch Proof Pin 1 mm Pin Touch Proof 1,5 mm DIN 42802
RoHS II	✓
REACH	✓