



ADHESIVE ELECTRODES FOR SLEEP OR LONG MONITORING



Description

Spes Medica adhesive electrodes have a special shape and a specific size studied to allow a perfect registration also on smaller muscles.

The innovative Hydrogel solid gel is developed especially for this application and gives to the electrode a stable adhesion on the skin. It also can be repositioned several times.

- Material used is suitable for PSG and highly durable.
- Small size for pediatric and facial use.
- Excellent recording performance due to gel application

UMDNS 11440

Materials

Electrode	AgAgCl
Material in contact with the skin	Hydrogel /Foam
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Ni with PE jacket
Liner	PET
Packaging	AL/PE

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	DENIFXXXXX
Dimensions	30x22 mm
Cable length	From 8 to 250 cm
Cable colour	○ ● ● ●
Connector	Touch Proof 1,5 mm DIN 42802

Environment

RoHS II	✓
REACH	✓