



BRIDGE ELECTRODE



Description

Spes Medica offers the gold standard in term of bridge electrode quality.

Our genuine bridge electrodes are worldwide famous for the quality of the recording, the durability of the materials, the perfect regulation and connection system.

The bridge is made by a plastic and hypoallergenic material which guarantees the maximum stability and the best patient's comfort.

The skin- electrode distance can be regulated in a safe and easy way thank to our special screw.

The screw is also a universal connector from plug, spring plug and alligator clip.

Our electrodes are available in adult and pediatric version and in 2 different materials: AgAgCl and sintered AgAgCl.

UMDNS 11440

Materials




Safety

Specifications

Environment

Electrode	AgAgCl
Bridge	ABS
Connector	Ag - Au
Electrode holder	PE
Packaging	PE

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Generic code	BEXXXXXXXXXX
Size	Adult, Pediatric, Neonatal
Bridge color	   Without bridge

PVC-Free	✓
RoHS II	✓
REACH	✓