



CHEST BAND



Description

Made by soft and elastic materials, it gives the best comfort to the patient. The band has easy snap-on connector for an easy connection with the cap tie. Measure indicator on the band.

UMDNS 10279

Materials

Band PA

Packaging PE

Safety

Biocompatibility
(UNI EN ISO 10993) Yes

Latex No

Phthalates No

Specifications

Generic code BELT000XXXX

Size & Colour code

From 40 to 50 cm (Pediatric)	●
70 cm (Adult)	●
95 cm (Adult)	●
120 cm (Adult)	●

Environment

RoHS II ✓

REACH ✓