




CONNECTION LEAD WIRES



Description

Our connection cables have been designed to minimize the noise during the connection to the electrode equipment.

The conductor is made by silver with Kevlar ensuring the maximum conduction and a high tensile strength.

The cable-connector connection is moulded to avoid the oxidation and any risk of breaking

UMDNS 16312

Materials

Lead wire	Cu/Ag with PVC jacket
Connector	CuZn/Au with PE jacket

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	BECXXXXXXXX
Patient side	Insulate Crocodile Mini Croco-Clip Spring 2-3-4 mm pin
Headbox side	1,5 mm Touch Proof DIN 42802 Spring 2 mm pin 2 mm pin
Cable length	50 - 250 cm
Cable color	

Environment

RoHS II	✓
REACH	✓