



# CUP ELECTRODE MRI CONDITIONAL



## Description

The cup electrode set is suitable for MRI tests ensuring safety conditions. For ICU the system has the advantage to avoid to remove the electrodes from patients' skin thanks to the easy disconnection of multipolar connector. These electrodes must be used with their dedicated connection cable.

UMDNS 17554

## Materials

Electrode	AgAgCl / Au
Material in contact with the skin	AgAgCl / Au
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn with PE jacket
Packaging	PP

## Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

## Specifications

Generic code	DCCCCPE-XX DGGG0GCE-XX
Electrode diameter	10 mm
Electrodes number	From 21 to 25
Recording Area	78 mm <sup>2</sup>

## Connection cable



Lead wire	Cu/Sn with PVC jacket
Connector	CuZn with PE jacket
Generic code	BECCHAR-XX
Connettore	Touch Proof 1,5 mm DIN 42802
Cable length	100 cm

## Environment

RoHS II	✓
REACH	✓