




# PEDIATRIC CUP ELECTRODE



## Description

The 6 mm diameter pediatric cup electrodes are made of thick and fine silver band avoiding bending. Thanks to the electrode dimension, the precision of recording is relevant on pediatric patients.

- Developed for pediatric patients and for recordings in more focused areas;
- They guarantee the same recording quality as the adult version;
- Pack of 10 colors.

UMDNS 17554

## Materials

Sensor	AgAgCl / Ag / Au
Lead wire	Cu/Ag with PVC jacket
Connector	CuZn/Au with PE jacket
Packaging	PP

## Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

## Specifications

Generic code	DXXXXXXXXX
Electrode diameter	6 mm
Recording Area	28 mm <sup>2</sup>
Cable length	From 100 to 250 cm
Cable color	
Connector	Touch Proof 1,5 mm DIN 42802

## Environment

RoHS II	✓
REACH	✓