



EAR CLIP ELECTRODE



Description

Ear clip electrodes are made with an over moulded spring clip for the best comfort of the patient as well as the best connection with the ear lobe.

- Usable as reference electrode;
- Available with several connectors and lead wire lengths.

UMDNS 11440

Materials

Electrode	Sn (on request AgAgCl, sintered AgAgCl or Gold)
-----------	---

Material in contact with the skin	PS
-----------------------------------	----

Lead wire	Cu/Ag with PVC jacket
-----------	-----------------------

Connection	CuZn/Au with PE jacket
------------	------------------------

Packaging	PP
-----------	----

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
--	-----


Latex	No
-------	----

Phthalates	No
------------	----

Specifications

Generic code	EARXXXXXX
--------------	-----------

Cable length	From 100 to 150 cm
--------------	--------------------

Cable color	
-------------	---

Connector	Touch Proof 1,5 mm DIN 42802
-----------	------------------------------

Environment

RoHS II	✓
---------	---

REACH	✓
-------	---