



# MONOPOLAR CABLES



## Description

Our connection cables have been designed to minimize the noise during the connection to the electrode equipment. The cable is made by silver with Kevlar ensuring the maximum conduction and a high tensile strength. The cable-connector connection is moulded to avoid the oxidation and any risk of breaking

UMDNS 16312

## Materials

Lead wire      Cu/ Ag with PVC jacket

Connector      CuZn/Au with PE jacket

## Safety

Biocompatibility  
(UNI EN ISO 10993)      Yes

Latex      No

Phthalates      No

## Specifications

Generic code      LEADXXXXXXXX  
MNCXXXXXXXX

Electrode side      Snap  
Spring Snap  
0,7 – 1 mm Pin  
1,5 mm Touch Proof

Headbox side      1,5 mm Touch Proof DIN 42802

Cable length      From 50 to 250 cm

Cable colour      ● ● ● ● ○ ●  
● ● ● ● ● ●

## Environment

RoHS II      ✓

REACH      ✓