




# SHIELDED MULTIPOLAR CABLES



## Description

Our connection cables have been designed to minimize the noise during the connection to the electrode equipment. The cable is made by silver with Kevlar ensuring the maximum conduction and a high tensile strength. The cable-connector connection is moulded to avoid the oxidation and any risk of breaking

UMDNS 16312

## Materials

Lead wire	Cu/ Ag with PVC jacket
Connector	CuZn/Au with PE jacket

## Safety

Biocompatibility <i>(UNI EN ISO 10993)</i>	Yes
Latex	No
Phthalates	No

## Specifications

Generic code	10XXXXXXX 20XXXXXXX
Electrode side	2-3 x 1,5 mm Touch Proof Pin 2-3 x Snap 2-3 x Pin 0,7 - 1 mm 2-3 x Insulated Crocodile
Headbox side	1,5 mm Touch Proof DIN 42802 5 pole DIN 45322
Cable length	From 100 to 300 cm

Cable colour	
--------------	--

## Environment

RoHS II	✓
REACH	✓