



SNAP ADHESIVE ELECTRODES



Description

Spes Medica adhesive electrodes have a special shape and a specific size studied to allow a perfect registration also on smaller muscles.

The innovative Hydrogel solid gel is developed especially for this application and gives to the electrode a stable adhesion on the skin. It also can be repositioned several times.

- High quality gel for the best adhesiveness;
- AgAgCl electrode for the best conduction;
- Easy snap connection

UMDNS 11440

Materials

Safety

Specifications

Environment

Electrode	AgAgCl
Material in contact with the skin	Hydrogel
Liner	PET
Packaging	AL/PE
Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No
Generic code	CDESXXXXXX
Dimensions	24 mm Ø 35x45 mm
Cable	no
Connector	Snap
RoHS II	✓
REACH	✓