



Électrodes Aiguilles / IOM




ANGLE 90°

TFDN351340	13 mm	● 0, 35 mm	100 cm	2,5 mm	10 pcs
TFDN351341	13 mm	● 0, 35 mm	150 cm	2,5 mm	10 pcs
TFDN351343	13 mm	● 0, 35 mm	250 cm	2,5 mm	10 pcs
TFDN451840	18 mm	○ 0, 45 mm	100 cm	5 mm	10 pcs
TFDN451841	18 mm	○ 0, 45 mm	150 cm	5 mm	10 pcs
TFDN451843	18 mm	○ 0, 45 mm	250 cm	5 mm	10 pcs
TFDN452540	25 mm	● 0, 45 mm	100 cm	5 mm	10 pcs
TFDN452541	25 mm	● 0, 45 mm	150 cm	5 mm	10 pcs
TFDN452543	25 mm	● 0, 45 mm	250 cm	5 mm	10 pcs
TFDN453740	37 mm	● 0, 45 mm	100 cm	10 mm	10 pcs
TFDN453741	37 mm	● 0, 45 mm	150 cm	10 mm	10 pcs
TFDN453743	37 mm	● 0, 45 mm	250 cm	10 mm	10 pcs
TFDN455040	50 mm	● 0, 45 mm	100 cm	10 mm	10 pcs
TFDN455041	50 mm	● 0, 45 mm	150 cm	10 mm	10 pcs
TFDN455043	50 mm	● 0, 45 mm	250 cm	10 mm	10 pcs

Particulièrement
adaptées aux
muscles faciaux



Positionnement des Aiguilles HOOK



Droit de l'abdomen

Hook
50 mm / 30°

Sous -cutanée
13 mm



Droit Fémoral

Sous -cutanée 13 mm

Hook 18 mm / 90°

Hook 25 mm / 90°

Hook 37 mm / 90°

