



Description

TechSystem is a Spes Medica innovative system to conduct EEG exams.

The kit is made by 3 different products, each one is necessary for EEG exam test:

- **TechPrep:** single use and single patient stick of degreasing gel for skin preparation.
- **TechDot:** single use and single patient conductive dots gel to be used as a conductive interface between the electrode and the skin to reduce the impedance.
- **TechTape:** single use and single patient pre-cut adhesive tape to be use for electrode fastening on the head of the patient.

TechSystem is a simple and innovative solution which allows to reduce:

- skin breakdown
- cross-contamination
- excess waste

The electrode is not included in the Tech System, any electrodes of your choice can be used

Composition

TechStick	Abrasive Gel
TechDot	Conductive Gel
TechTape	Adhesive Tape

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	KITECHXXX
N° TechStick	1 Stick
N° TechDot	24 Dots
N° TechTape	24 Tapes

Environment

RoHS II	✓
REACH	✓