



Description

The electrode for the pelvic floor can be used for both stimulation and recording.

It is an invasive device, previously placed on the examination glove, used by the operator on patients with pelvic floor diseases.

- Possibility to stimulate and record the pelvic floor;
- Biocompatible materials, latex free.

CONNECTION CABLE

Available different lengths and connectors for all the EMG equipment.

UMDNS 11441

Materials

Electrode	Ag
Support	PET or PU
Packaging	Medical paper/PP

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	PUDXXXX
Cable length	200 cm
Connector	Touch Proof 1,5 mm DIN 42802, 5 pole DIN45322

Connection cable

Materials

Cable	Cu/Sn with PVC jacket
Connector	PE
Packaging	PP

Product Specifications

Generic code	2001PUDXXXXXX
Cable length	200 cm
EMG Connector	5 poles DIN45322
STIM Connector	5 poles DIN45322; T.P. 1.5mm DIN42802



Environment

PVC-Free	✓
RoHS II	✓
REACH	✓