



## Description

The Urethral Electrode is made of an active part with an Au coating which allows a good quality in terms of stimulation and recording of the external urethral sphincter muscle during spinal, pelvic and urology surgery.

- Different catheter sizes available;
- 1 channel
- Touch proof connector 1.5 mm DIN 42802

UMDNS 11441

## Materials

Electrode	Ag with Au coating
Catheter	100% silicon
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Ni with PVC jacket
Packaging	Tyvek

## Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene Oxide (EtO)
Latex	No
Phthalates	No

## Specifications

Generic code	UEXXX
Number of poles	1 – 2
Catheter diameter	8-10-12-14-16-18 French
Catheter colour	<span style="color: brown;">●</span> <span style="color: black;">●</span> <span style="color: white;">○</span> <span style="color: green;">●</span> <span style="color: orange;">●</span> <span style="color: red;">●</span>
Cable length	200 cm
Connector	Touch proof 1,5 mm DIN 42802

## Environment

RoHS II	✓
REACH	✓