




BIPOLAR STIMULATION ELECTRODES



Description

Spes Medica has developed a moulded electrode with strain relief for saving the lead wire. Multipole shielded lead wire with high voltage rating resistance and strong handgrip for the best stability. Golded felts holder to avoid oxidation, new felts generation more stable with water. Double side Velcro band and a couple of stainless steel electrode and 10 replacement felts are included in the pack.

UMDNS 11441

Different configurations available on request.

Materials

Safety

Specifications

Environment

Sensor	Au
Material in contact with the skin	Stainless Steel / Felt
Handgrip	PE
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Au with PE jacket
Packaging	PP
Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No
Generic code	STILXXXX STISXXXX STIMXXXX
Electrode distance	10 mm (STIS) 25 mm (STIL) 24 mm (STIM)
Cable length	50-250 cm
Cable colour	Grey
Connector	Touch Proof 1,5 mm DIN 42802 5 pole DIN 45322
RoHS II	✓
REACH	✓

Accessories

Replacements felts



Replacements steel ball



Velcro fixing band



Generic code	FELTXXXX	SFERXXXX	VELCROBLUXX/3
Material	Felt	Steel	Velcro
Dimensions	Small - Large	Small - Large	From 30 to 100 cm