




# DISPOSABLE CUP ELECTRODE



## Description

The 10 mm diameter cup electrode is made by pure AgAgCl ensuring the best recording experience. The extra-flexible lead wire and the connector with strain relief guarantees the best signal transmission even for long - term monitoring.

- Highest record quality;
- Extra-flexible lead wire;
- Anti-break connector;
- 10 colors leadwire.

UMDNS 17554

## Materials

Electrode	AgAgCl
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Ni with PVC jacket
Packaging	AL/PE

## Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

## Specifications

Generic code	DAGDXXXXXX
Electrode diameter	10 mm
Recording Area	78 mm <sup>2</sup>
Cable length	From 100 to 250 cm Without cable
Cable color	
Connector	Touch Proof 1,5 mm DIN 42802 Pin 1 mm
Connection cable	DECSXXXXXX

## Environment

RoHS II	✓
REACH	✓