



DISPOSABLE HOOK NEEDLE ELECTRODES



Description

The Hook Needle Electrodes ensure excellent signal stability and optimal patient stimulation. They can be used when a long-time positioning of the needle in the same area is needed as during the Intraoperative Monitoring.

The hook guarantees a greater stability of the signal without risk of needle breakage, ensures a better contact with the skin and a lower impedance of the electrode contact.

UMDNS 11440







Materials

Needle	Stainless steel
Needle holder	PVDF
Needle cover	PVC
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Au with PVC jacket
Packaging	PP

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene oxide (EtO)
Latex	No
Phthalates	No

Specifications

Generic code	MNXXXXXXXXX
Needle diameter	0,35 - 0,40 mm
Needle length	12 - 40 mm
Cable length	From 8 to 250 cm
Cable colour	     
Connector	Touch Proof 1,5 mm DIN 42802 1 mm Pin

Environment

RoHS II	✓
REACH	✓