



GROUND BAND ELECTRODE



Description

The ground electrode ensures a high quality level of recording.

Spes Medica offers two different types of ground band: the elastic model ensuring the patient's best comfort and the pure steel version that can be used dry.

- 1 meter connection cable included in the pack

UMDNS 11441

Different configurations available on request

Materials

Sensor	Stainless Steel
Lead wire	Cu/Sn with PVC jacket
Connector	PE
Packaging	PP

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Nickel	No
Phthalates	No

Specifications

Generic code	GDRXXXXXXX
Size	From 28 to 60 cm
Connector	Touch Proof 1,5 mm DIN 42802

Environment

RoHS II	✓
REACH	✓