



DIGITAL RING ELECTRODES FOR RECORDING AND STIMULATION



Description

The pair of ring digital electrodes for stimulation and recording have cables marked with different colours (red and black) to indicate polarity.

They are available in:

- Version in stainless steel that ensures maximum durability and an area of conductive maximum precision;
- Paediatric and adult versions are available.

UMDNS 11441

Different configurations available on request.

Materials

Sensor	Stainless Steel
Electrode handgrip	PE
Lead wire	Cu/Sn with PVC jacket
Connector	Cu/Au with PVC jacket
Packaging	PP

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Nickel	No
Phthalates	No

Specifications

Generic code	SRRINGXXXXXX
Cable length	50-250 cm
Connector	Touch Proof 1,5 mm DIN 42802 5 pole DIN 45322
Size	Adult Pediatric

Environment

RoHS II	✓
REACH	✓