



Description

Spes Medica electrodes for central nervous system are disposable electrodes used for electrical stimulation and recording during the intraoperative monitoring.

The D-WAVE is a multipole array of electrodes that is positioned on the spinal column for the D-wave recording.

- Single use
- Distance between contacts optimized to get the best recording of the D-wave
- Flat recording area to get more stability during the recording
- Maximum flexibility

UMDNS 11453

Materials

Electrode	Platinum
Support	PU
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Ni with PVC jacket
Packaging	Medical Paper/PET

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene Oxide (EtO)
Latex	No
Phthalates	No

Specifications

Codes	IEP0184 IEP0185
Number of contacts	2 contacts (code IEP0184) 3 contacts (code IEP0185)
Contacts dimension	2x3 mm
Cable length	250 cm
Connector	Touch proof 1,5 mm DIN 42802

Environment

RoHS II	✓
REACH	✓