



Description

The electrode is conical and terminate in a perfectly blunted micro tip.

The isolated shaft is flexible and keep the shape for an easy modification by the operator. This probe is designed to reach small stimulation spots and can be adapted by the operator to reach any cavity.

- Light and ergonomic handgrip;
- Safety connector 1,5 mm DIN 42802.

UMDNS 11453

Materials

Electrode	Stainless Steel
Electrode insulation	PTFE
Handgrip	PS
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn with PVC jacket
Packaging	Medical paper/PET-PE

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene oxide (EtO)
Latex	No
Phthalates	No

Specifications

Generic code	SI1IFXXXXXXXXD
Tip diameter	0,5 mm
Probe length	95 mm
Cable length	250 cm
Connector	1,5 mm T.P. connector DIN 42802

Environment

RoHS II	✓
REACH	✓