



Description

SAC2 is an adhesive cream to be used with EEG Surface Electrodes. Provides resistance, cohesion and adhesiveness during use even if long term exams are made. SAC2 can be easily and quickly removed with water.

- Quick-drying cream
- No electrode damages
- Water soluble
- No residues
- No adverse smell
- External use only

UMDNS 11425

Composition

Ingredients

Water, Talc, Celite, Glycerol, CarboxyMethylCellulose, Sodium chloride, Potassium chloride, Phenoxyethanol, Ethylhexylglycerin

Tube material

Aluminum

Chemical and physical properties

pH

7,0 ÷ 8,7

Safety

Patch test on sensitive skin



Specifications

Generic code

SACX-XX

Packaging

10 tubes 100g each

Environment

RoHS II



REACH

